

SUPERA™ STENT IS LIKE NO OTHER STENT




**HIGH COMPRESSION
RESISTANCE¹**



**LOW CHRONIC
OUTWARD FORCE¹**



**HIGH FLEXIBILITY¹ AND
FRACTURE RESISTANCE²**



1. Test(s) performed by and data on file at Abbott.

2. Garcia, Lawrence A. et al., SUPERB Final 3-Year Outcomes Using Interwoven Nitinol Biomimetic Supera Stent. *Catheterization and Cardiovascular Interventions* 2017.